



When it's right, it just clicks

Philips Respironics Amara full-face mask - Performance, appeal, and simplicity at its best.

PHILIPS
RESPIRONICS
sense and simplicity

Amara - Philips Respironics full-face mask - performance, appeal, and simplicity at its best.



They say when things are right, they just click.

So the revolutionary **Amara** full-face mask from **Philips Respironics** is sure to click with medical professionals and patients alike. Its innovative, modular design is lighter, smaller, and has significantly fewer parts than leading traditional full-face masks. In fact, it requires just a single *click* to disassemble and reassemble the cushion and the mask frame for quick cleaning and replacement.

Philips Respironics has redefined the full-face mask, making it easy to fit and easy for patients to live with - plus it is designed and built to the exacting standards you would expect from Philips. **With Amara, it just clicks.**

We've raised our standard when it comes to performance, appeal, and simplicity.

For the medical professional...

Performance

- With only four sizes designed to fit over 95% of patient faces, Amara helps you get the right fit first time, every time.
- The flexible silicone cushion is designed to reduce pressure points, helping to keep your patients comfortable and compliant.

Appeal

- Smaller and lighter than other leading traditional full-face masks and designed with a high clarity finish - built to improve therapy acceptance.

Simplicity

- At just four parts, including an innovative forehead adjustment system, Amara makes light work of setting up your patients.
- With all cushion sizes being interchangeable between frames, Amara makes inventory management simpler and easier.

For the patient...

Performance


- As our best sealing and most comfortable traditional full-face mask yet, Amara is designed to enhance your quality of sleep.
- Customize and fine-tune the seal of your mask using the fine-glide forehead adjuster.

Appeal

- Amara was designed with you in mind. As the smallest leading traditional full-face mask on the market, it is less likely to get in the way of a good nights sleep.
- Less intimidating design.

Simplicity

- Amara's one-click cushion attachment gives you confidence that the mask is correctly assembled.
- With only four parts, assembly and cleaning are quick and easy so that you can spend more time doing the things you love.



One-click cushion assembly

A cushion design that reduces pressure points.

Fine-glide forehead adjuster

Eight selectable positions for fine-tuning seal and fit.

Lightweight minimal footprint frame

Smaller and lighter than any other leading traditional full-face mask on the market.

Quick-fit headgear clips

Reduces the need for continual headgear adjustment.

360° rotating elbow

Quiet-exhalation micro ports with quick release swivel.

Four-part modular design

Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.



Ordering information

| Description | Size | Part number |
|---|------|-------------|
| Amara mask with exhalation port | | |
| Amara mask with reduced size frame and reduced size headgear | P | 1090220 |
| Amara mask with reduced size frame and reduced size headgear | S | 1090221 |
| Amara mask with headgear | S | 1090222 |
| Amara mask with headgear | M | 1090223 |
| Amara mask with headgear | L | 1090224 |
| Amara mask without exhalation port | | |
| Amara mask SE with reduced size frame and reduced size headgear | P | 1090230 |
| Amara mask SE with reduced size frame and reduced size headgear | S | 1090231 |
| Amara mask SE with headgear | S | 1090232 |
| Amara mask SE with headgear | M | 1090233 |
| Amara mask SE with headgear | L | 1090234 |



For more information, visit the Amara web site. Simply scan the QR code with the QR reader app of your smart phone.

Contact us

Philips Healthcare is part of Royal Philips Electronics

How to reach us

www.philips.com/healthcare
healthcare@philips.com

Asia
+49 7031 463 2254

Europe, Middle East, Africa
+49 7031 463 2254

Latin America
+55 11 2125 0744

North America
+1 425 487 7000
800 285 5585 (toll free, US only)

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)

Philips Respironics International
Headquarters
+33 1 47 28 30 82

Philips Respironics Asia Pacific
+65 6882 5282

Philips Respironics Australia
+61 (2) 9947 0440
1300 766 488 (toll free, Australia only)

Philips Respironics China
+86 400 828 6665
+86 800 828 6665

Philips Respironics Deutschland
+49 8152 93 06 0

Philips Respironics France
+33 2 51 89 36 00

Philips Respironics Italy
+39 039 203 1

Philips Respironics Sweden
+46 8 120 45 900

Philips Respironics Switzerland
+41 6 27 45 17 50

Philips Respironics United Kingdom
+44 800 1300 845

www.philips.com/respironics

Respironics and Amara are trademarks of Respironics, Inc. and its affiliates.
All rights reserved.

Please visit www.philips.com/respironics



© 2012 Koninklijke Philips Electronics N.V.
All rights reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.
CAUTION: US federal law restricts these devices to sale by or on the order of a physician.
Printed in the United States

Broudy AE 8/3/12 MCI 4105039 PN 1102448