

Mirage Micro™ for Kids

NASAL MASK

The Mirage Micro™ for Kids is ResMed's second mask approved by the FDA and designed for pediatric patients¹ with sleep apnea. Based on the streamlined design and personalized comfort features of the original Mirage Micro, this mask combines headgear size and a child-sized, small

cushion to better fit children and their unique needs. The Mirage Micro for Kids is not only comfortable and easy for children to adjust to, but parents and sleep professionals will also find it easy to use.

MicroFit dial

Offers 24 positions for personalized fit, seal and comfort

Streamlined forehead support

Provides clear field of vision so therapy is less intimidating

Quiet, diffused vent design

Gently disperses air for uninterrupted sleep

360° elbow rotation


Multiple tubing positions for maximum comfort

Set-and-forget headgear clips

Easy for small hands to attach and detach; retain optimal headgear settings

Circular vent

Diffused air exits in 360° for improved comfort



Cushion loop

For easy removal of cushion



Mirage™ dual-wall cushion

Reduces pressure on nasal bridge and seals effectively; clipless for easy assembly and cleaning

Swivel

Easily attach and detach air tubing





TECHNOLOGY HIGHLIGHT

MicroFit dial

A choice of 24 positions at your fingertips, providing a comfortable, effective fit for children

- **Adjustable:** Can easily be adjusted with one hand while on the child's face
- **Precise:** Each turn of the dial adjusts the forehead support in small increments
- **Progressive:** Moves in toward the face, allowing the child to sense when the fit is most comfortable for them
- **Works in combination with the cushion:** As the dial adjusts the angle of the mask to the face, the dual-wall cushion adapts to enhance the seal without the need for overtightening the headgear

BENEFITS

Personalized fit and comfort

- MicroFit dial fine-tunes fit to suit the child's profile
- Mirage dual-wall cushion comes in a Small size, ideal for children 7 years and older

Quiet and streamlined

- Quietly and gently disperses air so child can experience uninterrupted sleep
- Slender design improves child's field of vision and reduces feelings of anxiety
- Detachable tubing is convenient if child needs to get up in the middle of the night

Intuitively easy to use

- MicroFit dial can easily be adjusted with one hand by the child or adult, even while lying down
- Fewer mask parts for easy cleaning and assembly
- Comprehensive user support package includes travel bag and user guide with mask fitting and care instructions for parents



- Child-friendly parts are easy for small hands to grasp and adjust



- Combination of MicroFit dial and dual-wall cushion achieves a perfect fit for a range of children's facial structures



- Mirage dual-wall cushion reduces pressure on small noses and sensitive nasal bridges while providing effective seal

PRODUCT CODES

Mirage Micro for Kids Nasal Mask

Complete System

61013